



**STAGNATION IS OUT
CHOICES ARE IN**

∨ **DISCOVER NOW**

WE LEARNED FROM THE BEST

Innovations based on feedback from patients and dedicated healthcare professionals like you.

As much as your patients love the fact that the Accu-Chek Insight system is intuitive and easy to use – we also received user feedback that the Diabetes Manager is not the fastest electronic device on the planet. We heard you and have been working hard on the next generation of improvements. We combined what our customers love most about our pump with new and exciting tools – giving you more reasons to recommend it to your patients. **So get ready for more choices!**

ROCHE DIABETES CARE OFFERS THREE UNIQUE WAYS FOR YOUR PATIENTS TO USE THE ACCU-CHEK INSIGHT INSULIN PUMP:

ACCU-CHEK® Insight
with **DBLG1®** from diabeloop



LOOP MODE

ACCU-CHEK® Insight
with mySugr®



APP-POWER

ACCU-CHEK® Insight



CLASSIC MODE

LOOP MODE

NOW AVAILABLE!

Accu-Chek Insight insulin pump
teams up with DBLG1 to offer
Automated Insulin Delivery (AID)
for your patients

DBLG1 handset
hosting the
algorithm

Discreet and reliable
Accu-Chek Insight
insulin pump

Dexcom G6 Continuous
Glucose Monitoring
(CGM) System

A SELF-LEARNING ALGORITHM BRINGS INSULIN DELIVERY TO THE NEXT LEVEL

ONE HANDSET HOUSES THE AID ALGORITHM AND OPERATES THE INSULIN PUMP

Based on real time glucose data, the algorithm does all the work in the background, deciding whether to adjust or even stop basal insulin delivery or deliver an automatic correction bolus – offering peace of mind for your patients.

THE ACCU-CHEK INSIGHT INSULIN PUMP

The first pump with a pre-filled NovoRapid® cartridge to be part of an AID solution.

THE DEXCOM G6 CGM

A small, wearable sensor and transmitter continuously checks your patient's glucose levels and delivers them to the smart DBLG1 algorithm every five minutes. With zero fingersticks and no calibrations*

QUICK START AND GO

Easy set-up will help your patients start using the hybrid closed-loop system in no time.

*If the glucose alerts and readings from the G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions.

WHAT IS THE DBLG1 SYSTEM?

- Three components – the Accu-Chek Insight insulin pump, the Dexcom G6 CGM and the DBLG1
- Diabeloop is an agile French MedTech company specialising in algorithm development and pioneering diabetes management



AUTOMATED INSULIN DELIVERY

Every five minutes, a glucose reading is sent to the handset via *Bluetooth®* technology. The algorithm analyses the data in real time, taking into account the patient's physiology, history and data entries (meals or exercise) to determine the correct dose of insulin to administer.

PERSONALISED & CONFIGURABLE MEETING INDIVIDUAL NEEDS

- Adjustable parameters that can be fine-tuned
- Automatic optimisation thanks to artificial intelligence

QUICK START WITH ONLY FOUR DATA POINTS

- Weight
- Total daily dose of insulin (TDD)
- Average amount of carbs per meal, in grams
- Safety basal rate profile

SAFE & SECURE

- Recommendation of rescue carbs and zen mode help avoid hypoglycemia
- When not in loop mode, a bolus can be given manually via the handset

EASY DATA MANAGEMENT

- AID data (CGM and pump) is stored on the web-based YourLoops platform – no cable or additional software needed
- Patients can easily grant access to their data

All your patient's diabetes data is collected in one place and automatically uploaded to the [cloud-based YourLoops platform](#).



LOOP MODE AT A GLANCE




AUTOMATED INSULIN DELIVERY

adjustment or even stopping
of basal rate, automatic
correction bolus



AUTOMATIC OPTIMISATION OF MEAL RATIOS

thanks to
auto-learning




PERSONAL ALERTS

and information
can be customised



6 CONFIGURABLE PARAMETERS

e.g. glycemic target,
low/high thresholds,
aggressiveness factors



USER INTERFACE

co-developed with
patients



SAFETY

Rescue carb recommendation
and zen mode ensure safety. If the
patient is not in loop mode,
a bolus can be given manually via
the handset.



Four clinical trials on adults* demonstrated the safety and performance of the DBLG1 System. Significant improvement in glycemic control was observed, along with a reduction of time in hypoglycemia while using the hybrid closed-loop system.

*Benhamou et al 2018: Benhamou PY, Hunecker E, Franc S, Doron M, Charpentier G; Diabeloop Consortium. Customization of home closed-loop insulin delivery in adult patients with type 1 diabetes, assisted with structured remote monitoring: the pilot WP7 Diabeloop study. Acta Diabetol. 2018;55(6):549-556.

APP-POWER

COMING SOON!

The Accu-Chek Insight insulin pump will be the first pump connected to the mySugr app



SOON YOUR PATIENTS CAN CONTROL THEIR INSULIN PUMP WITH ONE OF THE WORLD'S MOST POPULAR DIABETES APPS

HOW IT WORKS

Simply download the app from the Google Play store and set it up. Easy pairing of the pump and the mySugr app.* Open Bolus Advisor can be accessed via the mySugr app.

DISCREET BOLUS DELIVERY

Patients can discreetly deliver a bolus through the mySugr app on their Android smartphone.* All it requires is a personal authentication to ensure security.

ALL DATA IN ONE PLACE

Patients can manage their CGM, BG meter and pump data in one app.** All data gets transferred automatically in detailed reports.

SIMPLE DATA SHARING

Patients can print out mySugr reports or send them upfront via e-mail. No special software required. Quickly see the bigger picture.

*Compatible with Samsung Galaxy S7, S8, S9, S9+ and S10 smartphones

** For all mySugr app supported devices



APP-POWER AT A GLANCE



EASY PAIRING
of the pump and the
mySugr app



ALL DATA IN ONE PLACE
Patients can manage their
CGM, BGM and pump data
all in one app



DISCRETION
Give a bolus via the mySugr
app discreetly anywhere,
anytime



**OPEN BOLUS
ADVISOR**
directly accessible via
the mySugr app



FREEDOM OF CHOICE
Choose from two
pre-filled cartridges:
Fiasp® or NovoRapid®



SECURITY
Safe and secure insulin
delivery with personal
authentication



The use of the mySugr app has been shown to significantly lower the risk of severe hypoglycemic events, lower the estimated HbA1c and the risk of severe hyperglycemic excursions, improve blood sugar control, increase logging frequency by 10-20%, and trigger positive changes of glucose control.^{1,2,3,4}

CLASSIC MODE

If your patients are **IN** to a more minimal solution, then the classic pump is the right choice for them



SIMPLE, CLASSIC, EASY AND INTUITIVE

HOW IT WORKS

No need to carry additional devices, the pump can be used on its own – patients can adjust their insulin delivery directly, right on the pump. Zoom function, colour display and simple on-screen instructions for priming infusion sets and changing a cartridge make it easy.

DISCREET AND PERSONAL

The sleek and modern design is comfortable to wear and tucks away discreetly. Compatible with different infusion sets to meet individual needs.

SAFE AND RELIABLE

Offers various bolus options and the ability to run boluses in parallel. A very low basal rate suits individual insulin needs. Pause mode supports delivery of any missed insulin doses (after finishing a cartridge or changing a battery).

EASY AND CONVENIENT

Choose between two pre-filled cartridges, Fiasp or NovoRapid insulin – meaning there is no need for patients to fill by hand. Simple download of patient data with the Accu-Chek Smart Pix software for actionable insights and personalised care.



CLASSIC MODE AT A GLANCE




EASY AND INTUITIVE

Pre-filled insulin cartridges, zoom function and colour display



SAFE AND RELIABLE

Pause mode and various bolus options



SUPPORTING CHILDREN'S NEEDS

The classic pump option is particularly suited to children – with a very low basal rate and basal increments of 0.01 U/hr



CONNECTED

Data can be downloaded via the Smart Pix software



SELF POWERED

No need to carry additional devices



MODERN DESIGN

Comfortable and unobtrusive thanks to its sleek and modern design



EVEN MORE OPTIONS

New! Patients can choose between two pre-filled cartridges of Fiasp or NovoRapid insulin. Compatible with different infusion sets. Different cannula and tubing lengths to choose from.

WHO BENEFITS THE MOST FROM EACH CHOICE?



LOOP MODE

- Patients already using the Accu-Chek Insight insulin pump, who own or are considering a Dexcom G6 CGM, are interested in AID and are at least 18 years old
- Patients using other pumps who are due for contract renewal, already have CGM experience and are interested in AID
- Patients who would like to switch to an approved AID system

APP-POWER

- Patients (> 18 years) already using the Accu-Chek Insight insulin pump who want to control their pump more discreetly via their Android* smartphone
- Patients using other pumps who are due for contract renewal, and would like to control their pump remotely via the mySugr app
- Patients on MDI wanting to start insulin pump therapy and see the benefit of a pump controlled via an app on their own Android smartphone
- Patients who use other devices (CGM and BG meters)** and want to connect them all in one place

*Compatible with Samsung Galaxy S7, S8, S9, S9+ and S10 smartphones

**For all mySugr supported devices



THE CLASSIC

- Patients who are confident managing their diabetes therapy on their own
- Children who want an insulin pump might use the Accu-Chek Insight insulin pump as a stand alone option
- Patients who are not willing to use additional technology, no CGM, no AID
- Patients with good glucose control using an insulin pump and BGM/CGM





ACCUCHEK® Insight
with **DBLG1®** from  **diabeLoop**
LOOP MODE

ACCUCHEK® Insight
with  **mySugr®**
APP-POWER



ACCUCHEK® Insight
THE CLASSIC

1 Depong F, Mayer H, and Kober J. Real-World Assessments of mySugr Mobile Health App. *Diabetes Technology & Therapeutics*, June 2019. <http://doi.org/10.1089/dia.2019.0019>; **2** Hompesch M, Hergesheimer L, Kalcher K, Boubela R, Depong F. Retrospective Analysis of Impact on SMBG and Glycemic Control of Mobile Health (mHealth)-Application for Diabetes Management. Poster presented at DTM 2016, Bethesda, Maryland, USA; **3** Rose K.J, König M, Wiesbauer F. Evaluating success for behavioral change in diabetes via mhealth and gamification: mySugr's keys to retention and patient engagement. Poster presented at ATTD 2013, Paris, France; **4** Hompesch M, Scheiner G, Schuster L, Kober J, Depong F. Clinically relevant improvement in quality of blood sugar control in well controlled users of mySugr's mobile diabetes management tool. Poster presentation at the DTM 2018, Bethesda, Maryland, USA.

ACCUCHEK, ACCUCHEK INSIGHT, ACCUCHECK SMART PIX and MYSUGR are trademarks of Roche. NovoRapid®, Fiasp® and PumpCart® are registered trademarks of Novo Nordisk A/S. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by Roche is under license. DiabeLoop, DBLG1 and YourLoops are either trademarks or registered trademarks of DiabeLoop in Europe and other countries. Dexcom and Dexcom G6, are registered trademarks of Dexcom, Inc. in the United States and/or other countries. All other product names and trademarks are the property of their respective owners.

© 2020 Roche Diabetes Care

Roche Diabetes Care GmbH
Sandhofer Str. 116
D-68305 Mannheim